



INNOVATIONS THAT MATTER


---

## PRODUCT BROCHURE




In **Situ** Gas and Heating Systems

For **TEM** platforms



## Climate Nano-Reactor



The Climate Gas & Heating Solution enables dynamic studies of cyclic specimen transformations in gaseous environments at elevated temperatures and sub-Angstrom level.

Climate will convert your high vacuum TEM from a static imaging tool into a real-world research laboratory, enabling you to speed up your development of new catalysts or other energy relevant materials and techniques.

Climate is the only environmental solution in the market that allows full dynamic correlation of the structural and chemical data including reaction product analysis due to the integration with the optional dedicated DENSsolutions gas analyzer.


## Climate holder



## Typical applications



Heterogeneous catalysis



Nanomaterial growth &  
synthesis




Corrosion of  
metals & alloys



Green energy  
materials

## Selected publications




Time resolved in situ TEM images demonstrate the shape evolution of a Pd nanocrystals under 1 bar O<sub>2</sub> at 200 °C and the rebuilt Pd equilibrium structure by calculations. Scale bar: 5 nm.

### Shape evolutions of a metal nanocrystal in environmental conditions

The climate S3+ customers demonstrate an atomic scale TEM observation of shape evolutions of Pd nanocrystals under oxygen and hydrogen environment at atmospheric pressure. Combined with multi-scale structure reconstruction model calculations, the reshaping mechanism is fully understood.

These results give a direct insight into the behavioural response of nanoparticles to a 'real' reactive pressure environment, which is likely to improve the understanding of solid-gas reaction during catalytic applications.

Zhang, Xun, et al. "In situ TEM studies of the shape evolution of Pd nanocrystals under oxygen and hydrogen environments at atmospheric pressure." *Chemical Communications* 53.99 (2017): 13213-13216.




In Situ selected area electron diffraction (SAED) images. In Situ conditions: T = 600 °C, 200 mbar O<sub>2</sub>, 800 mbar N<sub>2</sub>, Flow = 0,6 ml/min..

### The reaction kinetics of nanoparticles

Acquiring the kinetics of gas–nanoparticle fast reactions under ambient pressure is a challenge owing to the lack of appropriate in-situ techniques. Now an approach has been developed that integrates time-resolved In Situ electron diffraction and a climate S3+ atmospheric gas cell system in TEM, allowing quantitative structural information to be obtained under ambient pressure with millisecond time resolution.

The ultrafast oxidation kinetics of Ni nanoparticles in oxygen was vividly obtained. This study gives new insights into Ni oxidation and paves the way to study the fast reaction kinetics of nanoparticles using ultrafast In Situ techniques.

Yu, Jian, et al. "Fast gas-solid reaction kinetics of nanoparticles unveiled by millisecond in-situ electron diffraction at ambient pressure." *Angewandte Chemie* (2018).



3D atomic structure reconstruction of Pt nanoparticles under the flow of selected gases.

In Situ conditions: T = 300°C, P = 1 Bar, Gas = O<sub>2</sub> (pure) or H<sub>2</sub> (5% diluted).


### Three-Dimensional Quantification of the Facet Evolution of Pt Nanoparticles in a Variable Gaseous Environment

The Climate G+ customers were able to quantify refaceting of Pt nanoparticles with atomic resolution during various oxidation–reduction cycles. In a H<sub>2</sub> environment, a more faceted surface morphology of the particles was observed with {100} and {111} planes being dominant. On the other hand, in O<sub>2</sub> the percentage of {100} and {111} facets decreased and a significant increase of higher order facets was found, resulting in a more rounded morphology.

This methodology opens up new opportunities toward in situ characterization of catalytic nanoparticles because for the first time it enables one to directly measure 3D morphology variations at the atomic scale in a specific gaseous reaction environment created by Climate G+.

Thomas Altantzis, et al. *EMAT, University of Antwerp, Belgium. Nano lett.* (2019).

## Why Climate?



### Experimental preparation made easy



#### 1. Full support with sample preparation

Experienced Application Engineers will provide you with tips and tricks on how to prepare your sample and a detailed workflow for FIB lamella preparation

#### 2. Predefine the experimental conditions

Calorimetry and Mass spectrometry data can be produced ex situ by the Climate system if required. This allows you to define the optimum experiment conditions before going to the TEM



### State of art environmental control

#### 1. Dynamic mixing

Using a specially designed and patented mixing valve instead of a premixing tank allows you to change the gas composition on the fly and to vary rates with an accuracy of 0.1 %

#### 2. Fast switching

The defined gas channel and minimal gas volume inside the Nano-Reactor enable you to change the gas environment within seconds

#### 3. Independent control


Gas composition, pressure and flow rate can be controlled independently using the widest experimental range

#### 4. Clean experiments


The modular design concept enables easy on-site exchange and cleaning of all critical components without gluing or welding

#### 5. Temperature accuracy and stability in gas environment

4-point probe heating provides the most accurate temperature with 0.01°C stability, even during gas flow. This allows accurate calorimetry data to be acquired during the experiment



### High impact results



#### 1. Data cross-correlation

Crystal structure, morphology, chemistry, thermodynamics and kinetics in one experiment

#### 2. Product analysis

The optimized gas analyzer and defined gas channel of the Nano-Reactor guarantees an accurate analysis of the environment at any time with a sensitivity of 5 ppm

#### 3. High Stability

Atomic resolution in TEM and STEM is routinely achievable in static and flow modes

#### 4. Optimized analytical capabilities

The optimized design enables EELS and large solid angle EDS collection

## Software for accurate environmental control

### Climate Software

Advanced control and data analysis

#### Guided system preparation

- A step by step procedure for the preparation of the system with multiple safety features included in the software, ensures a peace of mind

#### Stimulus control

- The 'Direct control mode' allows each user to set all relevant parameters for the experiment such as pressure, gas composition and the flow rate in the Nano-Reactor
- The 'Flowsheet' mode allows the advanced user to set every component of the gas supply system independently

#### Automate your experiments

- Design your experimental workflow and obtain reproducible results with an embedded profile builder


#### Safe operation

- The gas flow in the Nano-Reactor is monitored and includes an automatic safety check to stop the flow if required
- Flammability limit check
- Automatic data logging during the experiment



### Camera integration

True experimental integration



#### Easy & Fast experiments

- One SW to control in situ parameters & imaging

#### No mistakes or misinterpretation

- Synchronized in situ data with images at any frame rate

#### Save your time with your analysis

- Analyze images as a function of in situ parameters

## System specifications

JEOL	Climate Air	Climate G	Climate G+
Heating control	Four point probe resistive feedback		
Temperature range	RT to 1000 °C	RT to 1000 °C	RT to 1000 °C
Temperature Stability	0.01 °C	0.01 °C	0.01 °C
Pressure range	Ambient	0 - 1000 mbar	0 - 1000 mbar
Polepiece compatibility	All		
Alpha tilt range	UHR = ± 9 deg , FHP = ± 8 deg, HRP = ±15 deg, WGP = ±21 deg		
Modular design	Replaceable tubing and holder tip		
Resolution	≤100 pm*	≤100 pm*	≤100 pm*
Drift rate at in situ conditions	< 0.5 nm/min	< 0.5 nm/min	< 0.5 nm/min
Gas Mixing	N/A	Discrete	Continuous
Gas switching	N/A	< 60 s	< 1 s
Gas flow rate	Static	0 – 3 mln/min	0 – 3 mln/min
Micro-Calorimetry	V	V	V
EDS & EELS compatible	V	V	V

Thermo Fisher Scientific	Climate Air	Climate G	Climate G+
Heating control	Four point probe resistive feedback		
Temperature range	RT to 1000 °C	RT to 1000 °C	RT to 1000 °C
Temperature Stability	0.01 °C	0.01 °C	0.01 °C
Pressure range	Ambient	0 - 1000 mbar	0 - 1000 mbar
Polepiece compatibility	BioTWIN, S-TWIN, , TWIN, X-TWIN, C-TWIN		
Alpha tilt range	Bio-TWIN, C-TWIN, TWIN = ± 35 deg, X-TWIN, S-TWIN = ± 25 deg		
Modular design	Replaceable tubing and holder tip		
Resolution	≤100 pm*	≤100 pm*	≤100 pm*
Drift rate at in situ conditions	< 0.5 nm/min	< 0.5 nm/min	< 0.5 nm/min
Gas Mixing	N/A	Discrete	Continuous
Gas switching	N/A	< 60 s	< 1 s
Gas flow rate	Static	0 – 3 mln/min	0 – 3 mln/min
Micro-Calorimetry	V	V	V
EDS & EELS compatible	V	V	V

\*Listed specifications are dependent on microscope configuration

## Complete 'plug & play' package

1. Climate TEM specimen holder
2. Nano-Chips starter pack
3. Heating Control Unit
4. Laptop with pre-installed software
5. Chip alignment setup
- 6a. GSS light (Climate G)
- 6b. GSS (Climate G+)
7. Optional: Gas Analyzer

Including:  
Supporting tools



## Service and Support

Product Warranty	24 months with optional extension
Regulatory compliance	CE, RoHS, FCC
Mechanical compatibility	Approved by TEM manufacturers
Radiation safety	According to TEM manufacturers compliance regulations
Service	Dedicated Field Service Engineer




INNOVATIONS THAT MATTER

## DENSsolutions

Informaticaalaan 12  
2628ZD Delft  
The Netherlands  
Phone: +31 (0) 153 030 214  
Email: [info@DENSsolutions.com](mailto:info@DENSsolutions.com)  
[www.DENSsolutions.com](http://www.DENSsolutions.com)

To find your local distributor, please visit:  
[www.DENSsolutions.com/distributors](http://www.DENSsolutions.com/distributors)



Copyright © 2019 DENSsolutions B.V. All rights reserved.

Unless otherwise indicated, all materials on these pages are copyrighted by DENSsolutions B.V. All rights reserved. No part of these pages, either text or image may be used for any purpose other than personal use. Therefore, reproduction, modification, storage in a retrieval system or retransmission, in any form or by any means, electronic, mechanical or otherwise, for reasons other than personal use, is strictly prohibited without prior written permission.