For more, contact a Sales Engineer at sales@montanainstruments.com

CryoCore[®]

CryoAdvance™ 50

CryoAdvance™ 100

Cryostation® s200



| PERFORMANCE | | | | | NOTES |
|------------------------------------|-------------------------|--|--|--|--|
| Temperature Range | 4.9K-350K | 3.2K -350K | 3.4K-350K | 3.6K-350K | measured at platform |
| Sample Temperature Stability (P-P) | <100 mK | <10 mK | <15 mK | <20 mK | measured on manual positioner at base temp |
| Vibrational Stability (Platform) | <50 nm | <5 nm | <15 nm | <15 nm | peak to peak |
| Cool Down Time to 4.2K | ~3.5 hrs (to 6K) | ~2 hrs | ~3 hrs | ~10 hrs | typical performance |
| Cooling Power @ 4.2K | 100 mW (@ 6K) | 130 mW | 90 mW | 75mW | |
| OPTICAL PROPERTIES | | | | | |
| Optical Access | 5 optical ports | | | 9 optical ports | |
| | (4 radial + 1 overhead) | | | (8 radial + 1 overhead) | |
| Acceptance Angle | 60° full angle | 60° full angle | 27.4° full angle | 16° full angle | sample at center of platform |
| INTERFACING | | | | | |
| Electrical Access | 20 DC connections | 20 DC connections | 25 DC connections | 25 DC connections | |
| Interface Side Panels | Dual RF only | Dual RF + 25 DC (included) 1x configurable* | Quad RF (included) 5x configurable* | 6x configurable* | *RF, DC, fiber options available |
| DIMENSIONS | | | | | |
| Sample Space | ø53 mm x 63 mm | ø53 mm x 100 mm | ø100 mm x 116 mm | ø196 mm x 75 mm | |
| Beam Height | 89.9 mm | 140 mm | 140 mm | 100 mm | from table surface |
| OPTIONS | | | | | |
| Positioning | Manually adjustable | Manually adjustable, nanopositioner | Manually adjustable, nanopositioner | Manually adjustable, nanopositioner | |
| High NA Integration | None | Vacuum compatible objective vertically mounted in sample space | Vacuum compatible objective horizontally mounted in sample space | Vacuum compatible objective horizontally mounted in sample space | |
| Magnetic Field | None | 0.7T bipolar magnet integrated into a housing with 3 optical access ports. | Custom solutions required | Custom solutions required | |